AKUS™ ANAESTHESIA NEEDLE FOR SINGLE USE
Echogenic stimulation needle for nerve blocks. Two connection positions, Quincke 30°, detachable electrical cable with luer lock catheter.

INSTRUCTIONS FOR USE FOR NERVE BLOCK

Indications:

Needle with bevelled 25° tip designed for plexus and peripheral nerve blocks, using electrical nerve stimulation and/or ultrasound. Providing anaesthesia and analgesia in surgical procedures involving the upper and lower extremities, breast and abdomen. Useful needle for nerve block techniques in acute and chronic pain procedures.

The use of these types of needles is restricted to medical personnel with the training and appropriate experience necessary to perform nerve block techniques.

Instructions for use:

- Inform the patient about the technique being used and obtain written consent before starting the procedure.
- Closely inspect the proposed puncture site for infection and lesions before performing the puncture.
- Position the patient appropriately, according to the block being done. Correctly disinfect the area and cover with a sterile drape before starting the procedure. Inject a skin wheel of local anaesthetic at the puncture site.
- Inspect the puncture material (dosy), ensuring that it has been packaged correctly and the needle is sterile. *If not in perfect condition or sterility has been compromised, discard immediately.*
- Connect the transparent cable to the luer-lock connection on the distal end of the needle's body, purge the system with 0.9% Saline solution, or with the appropriate anaesthetic solution, until the liquid reaches the needle's tip.
- If performing the technique with nerve stimulation, connect the metal connection of the electrical cable to one of the two orifices found on the body of the needle. This can be done either on the left or right-hand side depending on the patient's position and/or the anaesthetologist's preference, and allowing connection of the nerve stimulator cable.
- Connect the other end of the electrical cable to the nerve stimulator and then the nerve stimulator cable to the skin electrode and place it on the body in an area near to the nerve being blocked. Set the initial current at 0.5 To 2 ma with a pulse width of 0.1 To 0.3 Ms.
- Punch the skin and slowly advance the needle towards the target nerve. If clear, muscular contractions are seen in the target nerve's area of distribution, reduce the electric current until the contractions disappear. If the level at which they disappear is between 0.3 And 0.5 Ma the stimulation needle is situated at an appropriate distance from the nerve. If muscular contractions are produced with a current less than 0.2 Retract the stimulation needle slightly as the tip is very close to, or within, the nerve fasciculate and this could cause lesions.
- If the technique is being performed under ultrasound, the akus™ needle features a system of specific grooves, in its last distal centimetre, allowing optimal, rapid visualisation of the needle tip. Likewise, the ergonomic grip ™ allows smooth, circumferential sliding movements making the ultrasonic tip flash more effective. This allows the technique to proceed safely and effectively as there is continual control of the tip position relative to the nerve.
- Puncture the skin and advance the needle slowly, using the ultrasound equipment to localise the desired nerve or nerve trunk.
- In both techniques, once correctly positioned and before administering the local anaesthetic, aspirate the syringe to ensure you are not close to a blood vessel.
- The akus™ plece needle's body is transparent which aids very rapid detection of blood. If this occurs, withdraw the needle and reposition it.
- Administer the desired anaesthetic, aspirating after every 5ml of anaesthetic solution.
- If desired the two techniques of nerve localisation can be combined: ultrasound and stimulation.
- The akus™ nerve block needle has length markings every 1cm. The ultrasound groove is only located in the distal centimetre.
- Withdraw the needle, disinfect the skin and place a sterile dressing.

Contraindications:

- Patient's refusal of the technique.
- Serious coagulation disorders, anti-coagulant treatment,
- Skin infection in the place of puncture.
- Septicaemia
- Allergy to local anaesthetics.
- Severe intracranial hypertension.
- Mild coagulations disorders.
- Known hypersensitivity to any of the materials used.
- Anatomic anomalies that make needle placement difficult via nerve stimulation and ultrasound.
- For patients with pacemakers or internal defibrillators, it is recommended to use ultrasound rather than electrical stimulation for nerve localisation.

Procedure complications:

- Local anaesthetics related toxicity by inadvertent intravascular injection, with neurological and /or cardiac symptoms.
- Nerve lesions caused by intrafascicular injection of local anaesthetic.
- Haematoma at the site of puncture.
- Pneumothorax
- Complete, extensive block due to inadvertent puncture of the epidural or intrathecal space, resulting in serious respiratory depression.
- Partial or incomplete block with insufficient anaesthetic due to poor technique.
- Infection (sepsis)
- Cardiac arrest
- Death

Warning:

- All akus™ nerve block needles have been sterilized by ethylene oxide.
- Single use needles and cables, should not be re-used.
- The akus™ needles is phthalate-free.
- Do not use if the packaging has been damaged.
- Ensure that the necessary material for general anaesthesia and cardiopulmonary resuscitation is readily available.

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