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**AKUSTM EPIDURAL KIT FOR SINGLE USE**
**Epidural anaesthesia kit: Tuohy epidural anaesthesia needle tip, 900 mm x 0.85 mm soft tip epidural catheter and adapter, connector, LOR loss of resistance syringe, 0.2µ filter, adhesive fixation strap for catheter, AKUS safety anti-puncture device, yellow epidural placement label, intrathecal catheter alert label.**

### INSTRUCTIONS FOR USE FOR THE EPIDURAL KIT

### Indications:

Localisation and placement of a catheter in the epidural space to perform epidural anaesthesia and/or analgesia.
The use of these types of needles is restricted to medical personnel with the training and appropriate experience necessary to perform epidural anaesthesia and/or analgesia techniques.

### Instrutions for use:

- Inform the patient about the technique being used and obtain written consent before starting the procedure.
- Inspect the proposed puncture site for infection and lesions before starting.
- Tattoos are present in this area, assess and explain the risks to the patient.
- Position the patient appropriately, according to the technique being done. Correctly disinfect the area and cover with a sterile drape before starting the procedure. Inject a skin wheel of local anaesthetic at the puncture site.
- Inspect the puncture material closely, ensuring that it has been packaged correctly and the needle is sterile. If not in perfect condition or sterility has been compromised, discard immediately.
- The epidural needle should be inserted with the stylet already inside it. It can be used with or without wings. The large and ergonomic AKUS™ wings allow the metallic cannula to be inserted with greater precision. Introduce 4-2-3 cm, using the cannula's centimetre surface markings for guidance, withdraw the stylet and attach the low resistance syringe, containing 3-4 centimetres of ether air or 0.9% saline solution. "Loss of resistance technique".
- Introduce the needle slowly until the epidural space is detected, using the loss of resistance technique with air or saline.
- Once the needle is ideally positioned, thread the epidural catheter through it. To aid this process you should use the small, insertion connector mounted on the catheter. For easier placement, remove the plastic pieces holding the catheter.
- Once the needle is in place in the correct position, insert the proximal end of the catheter into the click connector, this connector is showed openned. Insert the catheter up to the mark indicated on the proximal tip of the catheter, and will notice the limit, then close it and one click will indicate that the connector is perfectly sealed. The yellow connector is closed by pressing the upper cover on the lower base until the click is heard. The orange / blue connector closes with a 45° turn, until the seal is closed Both click connectors can be re-opened in case you want to remove the catheter for cut or reposition. The epidural kit is only supplied with one type of connector.
- After attaching the adaptor, do an aspiration test to ensure the catheter is not in the intrathecal space or within a blood vessel.
- It is recommended to use a bacterial filter. Inject 2mls of 0.9% saline solution on local anaesthetic to test the catheter's permeability.
- Disinfect the skin, place a sterile dressing and fix the catheter to the patient's back.
- Introduce the local anaesthetic solution through the catheter, using a bolus technique or continuous pump infusion.

### Kit contents:

- Tuohy epidural needle
- Catheter with adaptor
- Adhesive fixation strap
- Alus connector
- 0.2µ flat bacterial filter
- Low resistance syringe for detection of epidural space, LOR.
- Yellow sticker for correct signage of epidural technique.
- Red sticker warning of intrathecal catheter placement, if so desired.
- Alus safety anti-puncture device.

### Contraindications:

- Absolute:
  - Patient's refusal of the technique.
  - Serious coagulation problems, anti-coagulant treatment,
  - Skin infection in the place of punctation.
  - Septicemia
  - Allergy to local anaesthetics.
  - Severe intracranial hypertension.
- Relative:
  - Anatomic anomalies.
  - Skin tattoos.
  - Skin infection away from the proposed puncture site.
  - Mild coagulations disorders.
  - Non-corrected blood volume deficit.
  - Dehydration.

### Procedure complications:

- Inadvertent puncture of the intrathecal space.
  - PDPH (Post Dural Punction Headache) – post lumbar puncture headache
  - Hypotension
  - Urinary retention
  - Lumbago
  - Parasthesia
  - Partial paralysis
  - Total paralysis
  - Sudural or epidural haematoma
  - Meningitis
  - Systemic toxic reaction
  - Extensive block
  - Infection (epidural abscess)
  - Respiratory depression
  - Cardiac arrest
  - Death

### Warning:

- All AKUS™ epidural kits have been sterilized by ethylene oxide.
- Single use needle and material, should not be re-used.
- The AKUS™ needle is phthalate-free.
- Do not use if the packaging has been damaged.
- Ensure that the necessary material for general anaesthesia and cardiopulmonary resuscitation is readily available.

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**Zogear Industries Co., Ltd.**