**AKUS™ EPIDURAL KIT FOR SINGLE USE**

**epidural anaesthesia set: Tuohy needle, blunt ended epidural catheter 900mm x 0.85mm, low resistance syringe, 0.2µ filter, connector, adaptor.**

## INSTRUCTIONS FOR USE FOR THE EPIDURAL KIT

### Indications:

Localisation and placement of a catheter in the epidural space to perform epidural anaesthesia and/or analgesia.

*The use of these types of needles is restricted to medical personnel with the training and appropriate experience necessary to perform epidural anaesthesia and/or analgesia techniques.*

### Instructions for use:

- Inform the patient about the technique being used and obtain written consent before starting the procedure.
- Closely inspect the proposed puncture site for infection and lesions before starting.
- If tattoos are present in this area, assess and explain the risks to the patient.
- Position the patient appropriately, according to the technique being done. Correctly disinfect the area and cover with a sterile drape before starting the procedure. Inject a skin wheel of local anaesthetic at the puncture site.
- Inspect the puncture material closely, ensuring that it has been packaged correctly and the needles is sterile. If not in perfect condition or sterility has been compromised, discard immediately.
- The epidural needle should be inserted with the stylet already inside it. It can be used with or without wings. The large and ergonomic AKUS™ wings allow the metallic cannula to be inserted with greater precision. Introduce it 2-3 cm, using the cannula’s centimetre surface markings for guidance, withdraw the stylet and attach the low resistance syringe, containing 3-4 centimetres of either air or 0.9% saline solution. “Loss of resistance technique”.
- Introduce the needle slowly until the epidural space is detected, using the loss of resistance technique with air or saline.
- Once the needle is ideally positioned, thread the epidural catheter through it. To aid this process you should use the small, insertion connector mounted on the catheter. For easier placement, remove the plastic pieces holding the catheter.
- The soft silicone at the end of the catheter determines correct catheter insertion. This soft tip is designed to avoid perforation of the dura mater. The cannula has 3 soft boies positioned circumferentially around the catheter and has a radiopaque band along its entire length.
- The catheter’s position in the epidural space can be gauged from its surface markings. The first marking is placed at 5 cm, the double markings at 10 cm, the triple at 15 cm, and a quadruple at 20 cm. The distance between the markings is 1 centimetre.
- Once the catheter is ideally situated, withdraw the needle carefully ensuring the catheter’s position remains unchanged. It is expressly forbidden to remove the catheter through the steel cannula as the sharp bevel could damage and cut the catheter.
- On needle withdrawal, continue to control the catheter by monitoring these same markings.
- Once the catheter is ideally seated, insert its end into the lumen of the unscrewed Tuohy-Borst’s adaptor. Then twist until it secures the catheter firmly.
- After attaching the adaptor, do an aspiration test to ensure the catheter is not in the intrathecal space or within a blood vessel.
- It is recommended to use a bacterial filter. Inject 2mls of 0.9% saline solution or local anaesthetic to test the catheter’s permeability.
- Disinfect the skin, place a sterile dressing and fix the catheter to the patient’s back.
- Introduce the local anaesthetic solution through the catheter, using a bolus technique or continuous pump infusion.

#### Kit contents:

- Tuohy epidural needle
- Wings
- Catheter with connector
- Adaptor
- Flat bacterial filter
- Low resistance syringe for detection of epidural space, LOR
- Yellow sticker for correct signage of epidural technique.
- Red sticker warning of intrathecal catheter placement, if so desired.

### Contraindications:

- Absolute:
  - Patient’s refusal of the technique.
  - Serious coagulation problems, anti-coagulant treatment,
  - Skin infection in the place of puncture.
  - Septicaemia.
  - Allergy to local anaesthetics.
  - Severe intracranial hypertension.
- Relative:
  - Anatomic anomalies.
  - Skin tattoos.
  - Skin infection away from the proposed puncture site.
  - Mild coagulations disorders.
  - Non-corrected blood volume deficit.
  - Dehydration.

### Procedure complications:

- Inadvertent puncture of the intrathecal space.
- PDPH (Post Dural Punction Headache) – post lumbar punction headache
- Hypotension
- Urinary retention
- Lumbago
- Paraesthesia
- Partial paralysis
- Total paralysis
- Sudral or epidural haematoma
- Meningitis
- Systemic toxic reaction
- Extensive block
- Infection (epidural abscess)
- Respiratory depression
- Cardiac arrest
- Death

### Warning:

- All AKUS™ epidural kits have been sterilized by ethylene oxide.
- Single use needle and material; should not be re-used.
- The AKUS™ needle is phthalate-free.
- Do not use if the packaging has been damaged.
- Ensure that the necessary material for general anaesthesia and cardiopulmonary resuscitation is readily available.

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