KIT DE ANESTESIA EPIDURAL-ESPINAL (CSE) SOMENTE PARA USO ÚNICO

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SINGLE-SHOT COMBINED SPINAL-EPIDURAL (CSE) ANESTHESIA KIT

Epidural needle Tuohy 18G (80 mm), spinal needle Whitacre 27G (126 mm), soft tip epidural catheter 900mm x 0.85mm, loss-of-resistance syringe, 0.2µ filter, connector, adapter, transparent dressing for catheter fixation, yellow label for epidural access identification, red label for intrathecal access.

Indications:

Locate the epidural space, then carry out a spinal anaesthesia through the epidural needle, and finally place a catheter in the epidural space. A spinal- epidural technique is performed using a single lumbar space. Anaesthesia, and subsequent analgesia, are obtained in surgical procedures that require so. Anaesthesia for caesarean section with subsequent control of epidural analgesia, knee arthroplasty, hip arthroplasty, knee ligament-plasties, severe fractures of lower limbs (MMII) and pelvis, adjuvant abdominal surgery with subsequent control of analgesia, nephrectomy, radical prostatectomy, urological oncology surgery, abdominoplasty, vascular surgery for vascular endoprotheses.

Instructions for use:

- Before performing the technique, inform the patient and ask for their written consent.
- Before performing the puncture, rule out infections and injuries in the body area where the technique is to be performed.
- In case of tattooing of the area: Assess and explain the risks to the patient.
- Carry out a proper disinfection of the puncture site, and set up a sterile field prior to performing the technique.
- Observe the material carefully, ensuring the correct packaging and sterility of the kit. **In case of damage or loss of sterility, discard immediately.**
- Once the lumbar space where the puncture is to be performed has been chosen (it is recommended that the puncture be performed between the lumbar spaces L1 to L5), a small amount of local anaesthetic is administered at the puncture site.
- The epidural needle must be inserted with the stylet placed inside. Use the wings to aid insertion. The large, ergonomic AKUS™ wings allow for greater precision in the insertion of the metal cannula. Once it is 2-3 cm in, guided by the centimetre markings on the cannula, remove the stylet and connect the loss-of-resistance syringe, which can be used with 3 or 4 centimetres of air or physiological serum 0.9%. "Loss of resistance technique."
- Insert the needle slowly until the epidural space is located by losing resistance with air or serum.
- Once the needle has been placed in its correct position, the spinal needle is inserted through the epidural needle. This insertion must be carried out until the anterior prolongation of the body of the spinal needle is connected inside with the posterior hole of the body of the epidural needle. The spinal needle is designed so that once the bodies of both needles are perfectly assembled, the tip of the spinal needle protrudes 0.6 mm from the tip of the epidural needle. This distance is sufficient to access with the spinal needle to the subarachnoid space, allowing a very precise and safe control, as access to the epidural space without damaging structures or invading the epidural backspace.
- Once the subarachnoid space is located, remove the stylet and wait for CSF reflux. The transparent hub of AKUS™ introduces a magnifying glass , which facilitates the immediate observation of CSF. It's the only needle that offers this device. The magnifying glass and the size of the needle body, facilitate the proper inspection of normal CSF. If the CSF is not transparent, and presents turbidity or blood staining, assess aborting the technique, or perform puncture in different epidural space.
- The AKUS™ needle allows the use of LUER and LUER-Lock cone syringes. Once the syringe has been adjusted, with the local anaesthetic, press gently and insert it into the intrathecal space.
- Once the local anaesthetic has been administered, remove the spinal needle and insert the epidural catheter with the help of the small connector that comes already mounted with the catheter. Adjust this connector to the body of the epidural needle to facilitate the insertion of the catheter in the desired position. The tip of the catheter is determined by the coloured silicone at one end. This soft tip is designed to avoid perforation and unwanted introduction into the dura mater. The catheter has 3 outlet holes at the distal end of the catheter, and incorporates a radiopaque line along its entire length.
- The position of the catheter in the epidural space is determined by following the markings printed on it. The first single mark is located at 5 cm, the double mark at 10 cm, the triple at 15 cm, and the quadruple at 20 cm. As of that last marking, the separation from one marking to the next is 1 cm.

- After the catheter has been placed in the desired position, remove the needle without changing the position of the catheter. It is strictly forbidden to remove the catheter through the steel cannula, as the sharp bevel can damage and cut the catheter.
- Once the needle is removed and the catheter is positioned, the proximal end is inserted into the Tuohy-Borst adapter. This connector is now in the open position, so simply insert the catheter into the smallest hole in the adapter up to the marking at the proximal end of the catheter. Once inserted, make a 60° turn. This way, the catheter will be perfectly coupled with the adapter.
- After connecting the adapter, perform an aspiration test in order to avoid placing the catheter inside the intrathecal space (CSF reflux) or inside a blood vessel (blood reflux). The 0.2µ bacterial filter should be placed, and 2-3 ml of physiological serum should be administered to verify catheter permeability.
- Disinfect the skin by placing the dressing that comes with the kit, which serves to fix the catheter onto the skin of the back.
- Infuse anaesthetic solution through the catheter, using the bolus or continuous infusion techniques.

Contraindication:

- Absolute:
 - The patient does not consent to the technique.
 - Severe coagulation disorders, treatment with anticoagulants,
 - Infection of the skin at the puncture site.
 - Septicaemia
 - Allergy to local anaesthetics.
 - Severe intracranial hypertension..
- Relative:
 - Anatomical malformation
 - Skin tattoos.
 - Skin infections distant from the puncture site.
 - Slight coagulation alterations.
 - Uncorrected blood volume deficit
 - Dehydration.

Complications of the technique:

- Inadvertent puncture of the dura mater with the epidural needle with passage into the intrathecal space.
- PDPH (Post Dura Punction Headache) - post lumbar puncture headache
- Hypotension
- Urinary retention
- Low back pain
- Paraesthesia
- Partial paralysis
- Total paralysis
- Subdural or epidural haematoma
- Meningitis
- Systemic toxic reaction
- Central blocking
- Infection (epidural abscess)
- Respiratory depression
- Cardiac arrest
- Death

Warnings:

- All AKUS™ spinal needles have been sterilized with ethylene oxide.
- Single-use needle, should not be reused once used.
- The AKUS™ needle does not contain phthalates.
- Do not use if the packaging is damaged.
- Have the necessary equipment for general anaesthesia and cardiopulmonary resuscitation available.